TurboFlow Base Molds

TurboFlow Base Molds enable tissue to be positioned and embedded in wax when using routine or biopsy cassettes.

These re-usable base molds are made from stainless steel and the internal rounded corners allow for easy removal of the formed wax block. They can be used on all standard size cassettes. Available in packs of 12.

ltem	Size	Qty
TBM0707	7mm x 7mm	BX/12
TBM1515	15mm x 15mm	BX/12
TBM2525	25mm x 25mm	BX/12
TBM2530	25mm x 30mm	BX/12
TBM2537	25mm x 37mm	BX/12

